

Wysong Ergonomic Insole™

IMPORTANT INSTRUCTIONS

As you will note when examining these special insoles, they are very three-dimensional with many hills and valleys. This is not by accident. The anatomy of the foot and the dynamics of movement have been carefully studied, tested and incorporated into this design. When properly installed, the Wysong Ergonomic Insoles™ will bring your foot to life, create unequalled comfort and even ease and facilitate walking and running. You will literally feel like ‘going’ once you properly put these in your shoes. How the insole is placed in your shoe is critical to giving you these full benefits. Please follow these instructions carefully.

1. Universal size fits women’s size 5 to men’s size 15.

2. Remove shoe insole. The Wysong insole is substantial and may crowd the foot too much if you leave the shoe insole in place.




3. Place one of the Wysong insoles on the floor and place your foot on it to see how it should fit. The ends of your toes should just overlap the toe bridges on the insole. If you rock the foot forward you should feel the toes gripping and ready to push off forward. There, now you know how it should feel in the shoe.



4. You will note that unlike any other insole, the Wysong insole neither extends to the toe nor the heel in your shoe. It is not “too small.”

It is supposed to be short on both ends



5. If the insole is too large, the sides may be trimmed as well as the heel. When trimming, follow the contours. Do not trim from the end toe area.

Trim as necessary

